

DEBUG RECLINING ALL-TERRAIN WHEELCHAIR



DEBUG RECLINING ALL-TERRAIN WHEELCHAIR



The Reclining All-Terrain Wheelchair is designed to safely provide wheelchair users with transportation across areas where a conventional wheelchair becomes unmanageable. This includes areas such as beaches, snow covered sidewalks, forest trails, etc. Debug is intended to provide decades of service life in the harshest of environments.

LONG LASTING #316L STAINLESS STEEL FRAME / BACKREST ARTICULATES FROM UPRIGHT TO FULLY RECLINED POSITION / REMOVABLE TELESCOPIC LEG REST / FITS THROUGH A STANDARD ADA DOORWAY / LOW MAINTENANCE DELRIN BUSHINGS / WheelEEZ® POLYURETHANE BALLOON WHEELS ROLL ACROSS SAND, SNOW AND GRAVEL / DUAL STITCHED SLING WITH TENERA 130 THREAD / PHIFERTEX PLUS SLING MATERIAL RE-SISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE AND IS MACHINE WASHABLE.

SPECIFICATIONS

SEAT SIZES	ENGLISH (INCH)	METRIC (CM)	
Widths	14, 18, 22	36, 46, 56	
Depth	17	43	
Seat-to-floor height	18 – 21	43 – 53	
Backrest height	18	46	
Telescope footrest <i>knee to heel</i>	15 – 22	38 – 56	
Elevating leg rest <i>knee to heel</i>	15 – 22	38 – 56	
WEIGHTS			ARMREST STYLES
Capacity	250 – 400 lbs	114 – 182 kg	Removable swing away armrest
Chair weight	40 – 65 lbs	18 – 30 kg	Adjustable Height or Padded <i>optional</i> "Contoured" or specialized <i>optional</i> Interchangeable desk option available

Manufactured in the USA,
Debug has a 10 year warranty on stainless steel frames and a 30 day money back guarantee.
© Copyright Deming Designs Inc., 2001 – 2016. All rights reserved.

ORDER ONLINE WWW.DEBUGMOBILITYPRODUCTS.COM