

All Terrain Walker Conversion Kit Comparison Chart



WZ1-ATW-CK-24

WZ1-ATW-CK-24-2

Product Name	All Terrain Walker Conversion Kit 24	All Terrain Walker Conversion Kit 24-2
Price	See Item	See Item
Wheel Type	Low-Pressure Beach Balloon Wheel	Low-Pressure Beach Balloon Wheel
Best for use on	Soft Sand & Challenging Terrain	Soft Sand & Challenging Terrain
Maximum Payload	72 kg (158 lbs)	72 kg (158 lbs)
Axle Stem Legs	25.4 mm x 25.5 cm (1 x 10")	25.4 mm x 25.5 cm (1 x 10")
Axle Material	Marine Grade Anodized Aluminum	Marine Grade Anodized Aluminum
Wheel Dimension	24 x 12.3 cm (9.4 x 4.8")	24 x 12.3 cm (9.4 x 4.8")
Maximum Distance Between Brackets	33 cm (13")	33 cm (13")